



Provide the highest level of FiO₂ and dependable EtCO₂ detection for your high risk patients undergoing endoscopic conscious sedation procedures.

MEDIUM CONCENTRATION

LP/M	FiO ₂
8L/min	≈ 77%
12L/min	≈ 87%

MED/HIGH CONCENTRATIONS

LP/M	FiO ₂
High 10-15	≈ 88%-95%
Medium 8-12	≈ 77%-87%


PEDIATRIC USE

Min Flow Rate 4LPM
> 2 years of age
> 28 Lbs



For Patients With

- Obesity BMI>30
- OSA-Sleep Apnea
- Heart Disease
- Vascular Disease
- Respiratory Compromise or Disease



PREVENT. DETECT. RECOVER
 Questions? 855-766-0202